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XSEAT® light

Dekubitusprophylaxe-Sitzkissen

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operation manual

XSEAT® light

Pressure ulcer prophylaxis seat cushion

Please read the operating instructions carefully before using the system and heed the **warnings!**

Product name	XSEAT® light
scope of delivery	Decubitus trophylaxis seat cushion with foam core and PU cover (detachable by zip), instruction manual
Indications / areas of application	<ul style="list-style-type: none">• Pressure ulcer prophylaxis
Contraindications	Not known
Permitted patient weight min. - Max.	20-120 kg
material	<ul style="list-style-type: none">• Core: PU foams (HR / Visco)• Cover: PU-coated PES carrier material
Dimensions	See label on the product
Product weight	All sizes approx. 0.5 kg
Aid number (as defined in §139 SGB V)	
Medical device	Class 1, MDR 2017/745
Technical safety / flame retardancy	Flame retardancy test according to DIN EN 597 I + II, BS 7175 sec 3 crib 5, BS 7177
Testing for toxicological harmlessness / biocompatibility	DIN EN ISO 10993-5, DIN EN ISO 10993-1, DIN EN ISO 10993-10
Warnings	Before the first application, the system must be prepared for the first use (see "Preparing the system for use").
guarantee	1 year on the functional properties of the foam core (pro-rata guarantee)
maintenance	No
Cleaning instructions (pillowcase)	<p>The cover can be disinfected by wiping with commercially available disinfectants made from isopropanol (e.g. Bacillol) or with other suitable DGHM-listed products (e.g. ECOLAB cleaning cloth Incidin OxyWipe S XL). Chemothermal cleaning of the cover is possible. The cover can be reprocessed by autoclaving according to the guidelines of the RKI.</p>
Cleaning instructions (core)	<p>The foam core without cover can be cleaned with household means. This can be achieved by soaking the contaminated area or spraying the entire foam core with a disinfectant (eg Bacillol or ECOLAB foam spray Incidin OxyFoam S consisting of hydrogen peroxide). Then the foam core must be air-dried. The foam core can be reprocessed by autoclaving according to the guidelines of the RKI.</p>
Prepare the system for use	<ol style="list-style-type: none">1. Remove the pillow from the box.2. Place the pillow on the seat of the wheelchair or seat, "THIS SIDE UP" printed on the front upwards.

RKI = Robert Koch Institute

DHGM = German Society for Hygiene and Microbiology

MPG = Medical Devices Act

Autoclaving = fractional vacuum process at 105 degrees C and 55 minutes

Manufacturer in the sense of the MPG and MDR 2017/745

Funke Medical AG • Ährenfeld 10 • 46348 Raesfeld • Tel. : + 49-2865-266 93 50 • Fax: + 49-2865-8133

info@funke-medical.de • www.funke-medical.de

Download: 02.07.2020 13:10

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Funke Medical AG • Ährenfeld 10 • 46348 Raesfeld • Tel.: +49-2865-266 93 50 • Fax: +49-2865-8133