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## HEEL PROTECT®

Fersenfreilagerungsschuh

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## HEEL PROTECT®

Heel-free shoe

Please read the operating instructions carefully before using the system and heed the **warnings!**

<b>Product name</b>	<b>HEEL PROTECT®</b>
<b>scope of delivery</b>	Heel-free shoe with foam core and PU cover (removable with a zip), operating instructions
<b>Indications / areas of application</b>	<ul style="list-style-type: none"> <li>• Pressure ulcer prophylaxis</li> <li>• Pressure ulcer therapy</li> <li>• Therapy of the diabetic foot</li> <li>• Spitz foot prophylaxis (Pes equinus)</li> </ul>
<b>Contraindications</b>	Not known
<b>Permitted patient weight min. - Max.</b>	Suitable for all patient weights.
<b>material</b>	<ul style="list-style-type: none"> <li>• Core: PU foams (HR / Gel)</li> <li>• Cover: Mesh<sup>3</sup>, PU-coated PES carrier material</li> </ul>
<b>Dimensions</b>	Suitable for all shoe sizes.
<b>Product weight</b>	0.3 kg
<b>Aid number (as defined in §139 SGB V)</b>	11.11.05.0036
<b>Medical device</b>	Class 1, MDR 2017/745 
<b>Technical safety / flame retardancy</b>	Flame retardancy test according to DIN EN 597 I + II, BS 7175 sec 3 crib 5, BS 7177
<b>Testing for toxicological harmlessness / biocompatibility</b>	DIN EN ISO 10993-5, DIN EN ISO 10993-1, DIN EN ISO 10993-10
<b>Warnings</b>	<ul style="list-style-type: none"> <li>• Before the first application, the system must be prepared for the first use (see "Preparing the system for use").</li> <li>• Not suitable for changing patients</li> <li>• The shoe is only to be used in the supine position, not in the side or prone position.</li> </ul>
<b>guarantee</b>	1 year on the functional properties of the foam core (pro-rata guarantee)
<b>maintenance</b>	No
<b>Cleaning instructions (pillowcase)</b>	<p>The cover can be disinfected by wiping with commercially available disinfectants made from isopropanol (e.g. Bacillol) or with other suitable DGHM-listed products (e.g. ECOLAB cleaning cloth Incidin OxyWipe S XL). Chemothermal cleaning of the cover is possible. The cover can be reprocessed by autoclaving according to the guidelines of the RKI.</p> 
<b>Cleaning instructions (core)</b>	<p>The foam core without cover can be cleaned with household means. This can be achieved by soaking the contaminated area or spraying the entire foam core with a disinfectant (eg Bacillol or ECOLAB foam spray Incidin OxyFoam S consisting of hydrogen peroxide). Then the foam core must be air-dried. The foam core can be reprocessed by autoclaving according to the guidelines of the RKI.</p>
<b>Prepare the system for use</b>	<ol style="list-style-type: none"> <li>1. Remove the heel release shoe from the box upon receipt.</li> <li>2. The shoe is only to be used in the supine position, not in the side or prone position.</li> <li>3. When putting on the heel release shoe, it must be ensured that the heel is located above the recess and therefore does not rest anywhere.</li> <li>4. Tighten the Velcro fastener to prevent the heel release shoe from slipping.</li> <li>5. Adjust the hook-and-loop fastener for pointed foot prophylaxis (pes equinus) so that the foot is positioned at a 90 ° angle.</li> </ol>

RKI = Robert Koch Institute  
 DHGM = German Society for Hygiene and Microbiology  
 MPG = Medical Devices Act  
 Autoclaving = fractional vacuum process at 105 degrees C and 55 minutes

 **Manufacturer in the sense of the MPG and MDR 2017/745**

Funke Medical AG • Ährenfeld 10 • 46348 Raesfeld • Tel. : + 49-2865-266 93 50 • Fax: + 49-2865-8133  
 info@funke-medical.de • www.funke-medical.de

Download: 02.07.2020 12:42

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Funke Medical AG • Ährenfeld 10 • 46348 Raesfeld • Tel.: +49-2865-266 93 50 • Fax: +49-2865-8133