
Ref					
250700-250800	70-94cm	Ø 20mm	Ø 17mm	120kg	0.40 kg
243603-243604-243623-243624	70-94cm	Ø 22mm	Ø 19mm	120kg	0.36 kg
102023-102024	na	Ø mm	na	100kg	0.4 kg

243623 243624 243603 243604 102023 102024 250700 250800

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EN- Please read these notes carefully before use. Make sure to follow the instructions and keep them in a safe place ! If this aid is to be used by other persons, they must be fully informed of these safety warnings before use.

1. Safety indications

- ♥ ALWAYS consult with your physician or supplier to determine proper adjustment and use of the device.
- ♥ Don't store outside and don't let under the sun (Beware: don't let the product in a car in case of high outside temperature !!).
- ♥ Don't use in extreme temperature (above 38°C or 100°F or below 0° or 32°F).
- ♥ Check carefully all parts of this item (screws, tubes, tips, grips...).
- ♥ Check regularly if all parts are tightened (clips, handle, rubber tip..).
- ♥ Hold firmly and completely the handle when using.
- ♥ In case of falling : let the product fall down on the side, so that you don't fall on it.
- ♥ Always move forward quietly and with short steps to ensure a good balance of your weight and to avoid losing your balance.
- ♥ Don't use on wet or icy surfaces. (snow, ice, water..)!!
- ♥ Don't use in stairs, severe risk of injuries!!
- ♥ Use only as a walking aid and in proper way : don't use as a lever or to move furniture.
- ♥ ALWAYS observe the weight limit on the labelling of your stick.
- ♥ No modification must be made on this device. Use only original parts to replace the broken or defect one and let this be done by the technical staff.
- ♥ After assembling the walking aid, make sure that it is securely locked in OPEN position and level to the ground before using.
- ♥ ALWAYS test to see that the walking aid and attachments are properly

and securely locked in place before using.

- ♥ DO NOT hang anything on the Walking aid.
- ♥ In case you didn't use your walking aid for a long time, let it checked by the technical staff, before using.
- ♥ Be careful when using: beware holes and deformations in the floor.
- ♥ The suction cap must not be too worn down. If the cap shows considerable wear, ask your supplier or nursing staff to replace it.
- ♥ Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

2. Intended use, indications and contraindications

Intended use:

Walking aids are used for two purposes: as part of a rehabilitation programme when the user is recovering from an injury or operation and as a long-term aid to mobility when the user has a permanent difficulty with walking.

These devices help in maintaining balance, aid in the prevention of injuries (particularly falls), and allow the person to attend to his/her errands independently. Mobility aids can be used indoor or outdoor, depending on the model.

Mobility aids are used to modify the weight redistribution: some of the weight carried through the legs when walking is transferred through the arms of the frame or stick as it is leant on for support.

Indications

- Loss of autonomy due to walking disorders
- Balance and/or posture disorders
- Walking problems
- Rehabilitation programme after an injury or an operation

Contraindications

Don't use in case of:

- Severe cognitive disorder
- Severe loss of balance
- Severe motor or perception disorders
- Weakness in the upper limbs

Only specialists can adjust the product and teach you how to use it. A wrong use or adjustment could lead to overstress of the device or of the patient. The device must be adjusted individually by the technical staff to suit your needs. Ask the technical staff how to use this product and make sure that you understood clearly their explanation. In case of doubt, don't hesitate to ask again.

References : 250700,250800

Height Adjustment : Take out the clip from the tubes using both thumbs, then adjust the interior pipe to the desired height and lock again with the clip. Before using the crutch, check if the metal pin of the clip is correctly inserted into the pipe. The sticks have 1 adjusting hole on the outside pipe and 5 holes on the inside pipe, what means 5 adjusting positions.

The stick can be adjusted every 2 cm.

References : 243603,243604, 243623,243624

Height Adjustment : Take out the clip from the tubes using both thumbs, then adjust the interior pipe to the desired height and lock again with the clip. Before using the crutch, check that the metal pin of the clip has gone through both pipes and exits from the outer pipe.

The sticks have 1 adjusting hole on the outside pipe and 5 holes on the inside pipe, it means 5 adjusting positions.

The stick can be adjusted every 2 cm.

References: 102023, 102024

Height Adjustment: The correct cane length is the key to safe use and better mobility.

1. Obtain measurements while wearing regular walking shoes.
2. Standing upright, allow arms to relax (with normal bend at the elbow) at your sides.
3. Have a second person measure the distance from your wrist joint down to the floor. This number is the right length of cane for you.
4. Remove the rubber tip
5. Saw the cane at the correct length
6. Put the rubber tip on the stick.

To walk safely with a cane on level surfaces:

1. Hold the cane in the hand on your "good" side so that it provides support to the opposite lower limb
2. Take a step with the "bad" leg and bring the cane forward at the same time. Move the cane and affected leg forward together.
3. Lean your weight through the arm holding the cane as needed
4. Always have the bad leg assume the first full weight-bearing step on level surfaces
5. The cane should be moved the distance of one average step forward with

each move. You should not feel that you are stretching to catch up to the cane or stepping ahead of it.

If you are using the cane for general mobility rather than an injury, hold the cane using your dominant hand and bear weight on this side of your body. If you are working with a physical therapist due to an injury, he or she may have a specific cane-walking plan different from this one.

3. Cleaning

The stick has to be cleaned regularly with a mild cleaning agent and a soft cloth. NEVER use OILY substances!! Slip hazard!! Ask your retailer or the technical staff to check the walking aid regularly. The rubber tips must be checked and changed regularly.

4. Materials

Tube in anodised aluminium or wooden shaft and polypropylen handle.

5. Disposal

Dispose of this product in accordance with your regional waste disposal regulations.

6. Technical specifications

Check the table at the top of the user's manual. Duration of use: 2 years. This product has a 24 months warranty, from the date of purchase, against any manufacturing defect. Standards: DIN-EN-12182: 1999, EN-ISO 11334-1, DIN-EN-ISO 11334-4, DIN-EN-ISO 10993-1. Biocompatible.

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