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Hyper Foam PLUS clinic

Dekubitustherapie-Matratze

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Hyper Foam PLUS clinic

Pressure ulcer therapy mattress

Please read the operating instructions carefully before using the system and heed the **warnings!**

Product name	Hyper Foam PLUS clinic
scope of delivery	Decubitus therapy mattress with foam core and PU cover (removable with a zip), operating instructions
Indications / areas of application	<ul style="list-style-type: none"> • Pressure ulcer prophylaxis • Pressure ulcer therapy up to grade 3 (Seiler / EPUAP) • Pain therapy • For patients with medium risk of pressure ulcers
Contraindications	Not known
Permitted patient weight min. - Max.	20-140 kg
material	<ul style="list-style-type: none"> • Core: PU foams • Cover: PU-coated PES carrier material
Dimensions	Width: 80-100 cm Length: 190-220 cm
Product weight	11 kg (90/200 cm)
Aid number (as defined in §139 SGB V)	
Medical device	Class 1, MDR 2017/745
Technical safety / flame retardancy	Flame retardancy test according to DIN EN 597 I + II, BS 7175 sec 3 crib 5, BS 7177
Testing for toxicological harmlessness / biocompatibility	DIN EN ISO 10993-5, DIN EN ISO 10993-1, DIN EN ISO 10993-10
Warnings	<ul style="list-style-type: none"> • Before the first application, the system must be prepared for the first use (see "Preparing the system for use"). • Suitable for adjustable lying surfaces • Side rails at least 22 cm above the top of the mattress
guarantee	3 years on the functional properties of the foam core (pro-rata guarantee)
maintenance	No
Cleaning instructions (mattress cover)	<p>The cover can be disinfected by wiping with commercially available disinfectants made from isopropanol (e.g. Bacillol) or with other suitable DGHM-listed products (e.g. ECOLAB cleaning cloth Incidin OxyWipe S XL). Chemothermal cleaning of the cover is possible. The cover can be reprocessed by autoclaving according to the guidelines of the RKI.</p>
Cleaning instructions (core)	<p>The foam core without cover can be cleaned with household means. This can be achieved by soaking the contaminated area or spraying the entire foam core with a disinfectant (eg Bacillol or ECOLAB foam spray Incidin OxyFoam S consisting of hydrogen peroxide). Then the foam core must be air-dried. The foam core can be reprocessed by autoclaving according to the guidelines of the RKI.</p>
Prepare the system for use	<ol style="list-style-type: none"> 1. Remove the mattress from the box immediately upon receipt. 2. Put the mattress in the bed. 3. Note the print on the top of the mattress "THIS SIDE UP".

RKI = Robert Koch Institute
 DHGM = German Society for Hygiene and Microbiology
 MPG = Medical Devices Act
 Autoclaving = fractional vacuum process at 105 degrees C and 55 minutes

Manufacturer in the sense of the MPG and MDR 2017/745

Funke Medical AG • Ährenfeld 10 • 46348 Raesfeld • Tel. : + 49-2865-266 93 50 • Fax: + 49-2865-8133
 info@funke-medical.de • www.funke-medical.de

Download: 02.07.2020 12:09

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Funke Medical AG • Ahrenfeld 10 • 46348 Raesfeld • Tel.: +49-2865-266 93 50 • Fax: +49-2865-8133