PA100, PB100, PC100, PD100, PF100

EN Instruction for Use

1. Safety information

⚠ WARNING: The probe is to be operated by trained personnel only.

⚠ WARNING: Do not use the probe and SpO2 device in the presence of flammable anesthetics or gas to prevent explosion hazard.

⚠ WARNING: Do not use the probe and SpO2 device in the Magnetic Resonance Imaging (MRI) ambience.

⚠ WARNING: The SpO2 readings and pulse signals can be affected by the conditions of ambience, probe, and patient.

⚠ WARNING: Do not disassembly the probe. The SpO2 device is without any user-serviceable part inside and only qualified service personnel can perform maintenance service.

⚠ WARNING: Do not expose the probe and SpO2 device to extreme moisture (such as rain) to ensure accurate performance and device safety.

⚠ WARNING: If the accuracy of measurement by the SpO2 device is uncertain, check the patient's vital signs by alternate means.

⚠ WARNING: This probe and SpO2 device is intended only as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms.

⚠ WARNING: Reposition the probe at least once every 4 hours to allow the patient's skin to breath and to check patient's condition regularly.

⚠ WARNING: Do not use the SpO2 device with other devices (such as, the cuff of blood pressure monitor) that may interfere with blood flow and cause inaccurate measurement.

⚠ WARNING: The SpO2 device will be affected by electromagnetic interference or strong ambient light source during operation.

⚠ WARNING: User should stay calm and position finger stably. The accuracy of measurement taken right after exercise or during hand shaking could be compromised.

⚠ WARNING: The probe and SpO2 device should not be applied to a body part other the finger or to a wounded body part

⚠ WARNING: Do not overextend the device's spring.

⚠ WARNING: Do not look directly inside the housing during the measurement. The red light and the invisible infra-red light in the probe are harmful to your eyes.

2. Introduction

The probe combined with pulse oximeter is to spot-check oxygen saturation in blood (SpO2) and pulse rate.


The different probes combined with pulse oximeter is used on adults, neonatal, pediatric at hospital, clinics, and/or home, and disposable and extended probe.



The probe contains a dual light source and a photo-detector.

3. Features

- Measure and display reliable SpO2 value and pulse rate.
- High-brightness LED design

4. Product introduction

PA100		PB100	
			
1. Handler	2. Probe	1. Probe	2. Connector
3. Connector	4. Cable	3. Cable	

PC100		PD100	
			
1. Handler	2. Probe	1. Probe	2. Connector
3. Connector	4. Cable	3. Cable	

PF100	
	
1. Probe	2. Connector
3. Cable	

5. Operation

Press handler to open up the probe. Put one finger in the middle of probe and release the handler.

If one finger can not fit-in, switch to another one.

You may pull up the probe to enlarge the space for various fingers.

Connect to SpO2 device to start measurement.

⚠ WARNING: Please confirm the type of SpO2 device. If it is not the type, that might cause damage.

⚠ WARNING: As with all medical equipment, carefully route patient cabling to reduce the possibility of patient entanglement or strangulation.

⚠ WARNING: Do not reuse or disinfect disposable SpO2 probe. (only for PD100).

⚠ WARNING: It may damage the efficiency of the device if not apply with a Rossmax compatible oximeter.

6. Specifications

6.1 Performance

Scope of measurement: SpO2: 35% - 100%

Pulse rate: 30-250 bpm (beats per minute)

Accuracy: SpO2: 70%-100% : ±2%,
35% - 69%: unspecified

Pulse rate: 30-250 ± 3 bpm

6.2 Electrical specifications

Red LED wavelength: 660/680 nm

infrared LED wavelength: 880/905/940 nm

Connector Type: DB 9-pin

6.3 Environmental conditions

Operating temperature 5°C – 40°C (41°F – 104°F)

Storage temperature: -25°C – 70°C (-13°F – 158°F)

Relative humidity: 15% -90% (no condensing)

6.4 Physical characteristics

PA100 (for Adult)

weight:64g, Size:68x26x30mm, Length:1000mm

PB100 (for Neonatal)

weight:34.4g, Size:13x11x5.6mm, Length:1000 mm

PC100 (for Pediatric)

weight:34.6g, Size:29.4x17x20mm, Length:1000 mm

PD100 (for Adult)

weight:16g, Size:13.5x80.5x2.0mm, Length:1000 mm

PF100 (for Adult)

weight:40g, Size:48x39x19mm, Length:1000 mm

6.5 Optional Accessory

Extension cable(PE100)

weight:92.8g, Size:45.5x29.7x21.4 mm, Length:2400 mm

6.6 Standards

IEC60601-1-2, Class B; IEC60601-1, Type BF;

ISO 80601-2-61

6.7 Markings

	ATTENTION
	CE Mark
	Dispose on your local recycling center for safe treatment.
	Do not re-use for single use probe
	Date of manufacture
	Manufacturer

	Use-by date
	CONSULT INSTRUCTIONS FOR USE

7. Problem shooting and maintenance

7.1 Cleaning

Surface cleanings is by using a soft cloth dampened with either a commercial, non-abrasive cleaner or a solution of 70% (w/w) Ethanol alcohol in water, and lightly wiping the surfaces of the oximeter probe. Clean the LED and photo-sensor with moist cloth or cotton ball and alcohol gently.

The aforementioned general cleaning process is not for infection prevention. Please contact the specialist for the process of contagious infection.

7.2 Disposal

At the end of its life, the probe should not be disposed of in household rubbish. Enquire about the options for environment-friendly and appropriate disposal. Take into account local regulations.

8. Electromagnetic interference

Caution: This probe has been tested and found to comply with the limits for medical devices to the IEC 60601-1-2 and MDD 93/42/EEC. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare environments (for example, electrosurgical units, cellular phones, mobile two-way radios, electrical appliances, and high-definition television), it is possible that high levels of such interference due to close proximity or strength of a source may result in disruption of performance of this device. This probe is not designed for use in environments in which the pulse can be obscured by electromagnetic interference. During such interference, measurements may seem inappropriate or the monitor may not seem to operate correctly.

9. Warranty

The company warrants pulse oximeter at the time of its original purchase and for the subsequence time period of one year. The company warrants SpO2 probe free of defects at the time of its original purchase and for the subsequence time period of six months.

The warranty does not cover the followings:

- The device series number label is torn off or cannot be recognized.
- Damage to the device resulting from misconnection with other devices.
- Damage to the device resulting from accidents.
- Changes performed by users without the prior written authorization of the company.

 Rossmax Innotek Corp.
12F., No. 189, Kang Chien Rd., Taipei, 114, Taiwan.
Tel: +886-2-2659-7888
Fax: +886-3-2659-7666
www.rossmax.com

The text is subject to change without further notice.

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