Pedaltrainer RFM® digital

DEUTSCH

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
ENGLISH

1. Proper use:

The **Pedal Trainer RFM® digital** is intended to be used to improve the coordination and circulation in the extremities. It serves to build up the musculature and train motion patterns.


2. Warning:



Read the operating instructions before you use the **Pedal Trainer RFM® digital** for the first time and store it. Make sure that the **Pedal Trainer RFM® digital** is in a stable position and use it regularly. Increase duration and intensity slowly and consult your physical therapist or doctor if you feel unwell or feel pain. **The maximum load is 50 kg!**

3. Set-up/ Assembly:

- 1.) Carefully remove the **Pedal Trainer RFM® digital** from the shipping box (see Fig. 1) and remove the red safety strip at the display.
- 2.) Unfold the support legs of the **Pedal Trainer s RFM® digital**, until they audibly click into place; set it on a stable, even surface and screw the adjustment screw A into place (see Fig. 2)
- 3.) Turn the adjustment screw A, in order to adjust the pressure for axis B. Turning to the left increases the pressure; turning to the right reduces it. Set it back to minimum when you are finished using it. Press switch C to activate the built-in computer. Pressing the button repeatedly pages through the various modes (see Fig. 3)

- 4.) Pull out the locking pin in order to unlock the **Pedal Trainer RFM® digital** after use to fold it up and stow it away, saving space. (see Fig. 4 and 5)
- 5.) The counter with display – for “reset” please press button C for a few seconds.

<p>Time: Press button C 1x : arrows next to CNT and Time</p> <p>The device displays the total rotations and the elapsed time. It does not stop when you stop pedalling but continues to run!</p>	<p>Calorie consumption: Press button C 2x: Arrows next to CNT and CAL</p> <p>The device displays the total rotations and die calorie consumption.</p>
<p>UPM: Button C 3x press : Arrows next to CNT and RPM</p> <p>The device displays the total rotations and die rotations/minute.</p>	<p>Scan Modus: Button C 4x press : Arrows next to CNT and RPM and Scan</p> <p>The device displays the total rotations and alternately the elapsed time, calorie consumption and rotations/minute. Switches after about 10 seconds.</p>

The device shuts down after around 3 minutes' idle time.

4. Safety warnings:

- (a) Prior to use check to make sure that the device is in proper condition, the supporting legs are locked in place and the Pedal Trainer is standing securely with all 4 feet on the ground.
- (b) Set the desired pressure again before each use.
- (c) Longer use can cause the bottom bracket to heat up.

5. Cleaning instructions:

Clean the **Pedal Trainer RFM® digital** with a damp cloth and at regular intervals (absolutely necessary when you change users) using a disinfectant. You can disinfect it with regular disinfectants commercially available following the manufacturer's instructions, e.g. RFM Sept. Do not use any abrasive or corrosive cleaners or high pressure or steam cleaners!

6. Changing the battery:

Carefully press the little clip on the side in using a small screwdriver, until the LCD-unit (display) comes loose. Remove the LCD-unit.

Behind it you will find the control unit with the button cell battery (LR44). Remove it from the control unit.


Using a small screwdriver, replace the battery (LR44) and set the control unit back in the housing. Place the LCD-unit in the front again and press it into the frame until it clicks into place.

Daimlerstr. 12 a
D-25337 Elmshorn
Tel.: (0)4121 4733 12
Fax: (0)4121 4733 23
www.rehaforum.com
info@rehaforum.com

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Information concerning the disposal of old electric and electronic products

The pictured symbol on the product, the packing or the enclosed documents signifies that the product should not be disposed of after the end of its life cycle together with the domestic waste. You as consumer are obliged by law to render electric and electronic products at a local place of collection for electronic equipment. This is at no charge. The proper disposal enables the recycling of valuable resources and avoids possible negative effects on humans and the environment. For further information please contact your local authority or the vendor of the product.

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