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## EN FLEXIBLE LIPIDO-COLLOID CONTACT LAYER

### PRESENTATION

**URGOTUL** incorporates Technology Lipido Colloid (TLC) exclusively developed by Urgo. **URGOTUL** is a non-occlusive flexible and conformable dressing comprising of a polyester mesh impregnated with hydrocolloid particles dispersed in a petroleum jelly matrix. This mass contains carboxymethyl cellulose (hydrocolloid), paraffin oil, petroleum jelly and carrier polymers.

### PROPERTIES

In contact with wound exudate, the hydrocolloid particles in the **URGOTUL** mass interact with the petroleum jelly component to form a lipido-colloid interface. The interface creates a moist environment which is ideal for wound healing.

Due to the chemical composition, **URGOTUL** is non greasy to the touch and has non-slip properties. **URGOTUL** will not adhere to either the wound or the surrounding skin. Removal of the dressing is painless and atraumatic for the patient. As a result, Urgotul offers some specific properties:

- Creating and maintaining a moist environment that promotes healing,
- Gels in contact with the wound,
- Removal which does not damage newly formed tissue,
- Pain-free dressing changes for the patient.

### INDICATIONS

**URGOTUL** is indicated for the management of:

- Acute wounds (burns, traumatic wounds, abrasions, post-operative wounds) and chronic wounds (leg ulcers, pressure ulcers and diabetic foot ulcers) at the granulation and epithelialisation stage.
- Epidermolysis bullosa (EB).
- Flexible and very conformable, Urgotul is especially suitable for difficult to dress or awkwardly placed wounds and for packing deep wounds.

### INSTRUCTIONS FOR USE

- Prepare the wound according to local protocol. If an antiseptic is first used, rinse the wound thoroughly with saline solution before applying **URGOTUL**.
- Remove the protective transparent film.
- Apply **URGOTUL** to the wound in a single layer.
- **URGOTUL** can be cut with sterile equipment to fit the size of the wound and the surrounding skin.
- Cover **URGOTUL** with a secondary dressing suitable to absorb the level of wound exudate.
- Hold secondary dressing in place with a conforming bandage (K-Band / K-Lite), adhesive tape or an elasticated tubular bandage.
- **URGOTUL** should be changed every 2 to 4 days according to the type of wound and may be left in place for up to 7 days dependent on wound clinical condition.
- For patients with epidermolysis bullosa, the dressing should be changed after 1 to 3 days.

### PRECAUTIONS FOR USE

- **URGOTUL** adheres to latex surgical gloves. It is therefore recommended that gloves be moistened with saline in order to make the handling of Urgotul easier.
- If there are clinical signs of local infection and depending on clinical judgement, the wound can be treated with an antibacterial dressing such as Urgotul SSD/Urgotul Silver//Urgotul Duo Silver.
- In the case of deep wounds and fistula wounds, a section of the dressing should be left visible to enable easy removal.
- Do not re-sterilise the dressing.
- Single use sterile individual packaging: re-using a single use dressing may lead to risks of infection.

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