

Infusion - L86 - vented sets

L86 (closed cap, MLL)



1. The special design of the spike enables easy penetration of container ports and offers a safe connection to the container. Vented version. Air vent with hydrophobic, sterile air filter and ventilation cap. Standard drip tube (20 drops/ml). The protection cap maintains sterility of the set.
2. The flexible, transparent non-reflective material of the drip chamber with level ring mark enables easy priming and control of drip rate. Equipped with 15 micron fluid filter.
3. Precision roller clamp enables exact adjustment of the desired flow rate maintaining it with minimum variation. Additionally, it is equipped with a clamp/hanging device for tubing and for protection against stick accidents after use of the set.
4. Transparent, flexible PVC tubing with plasticiser TOTM 3,0/4,1 mm (ID/OD) in the set (regular bore).
5. Male Luer-Lock adapter guarantees a safe connection to the venipuncture device. The FlowStop protection cap for Male Luer-Lock adapters with integrated hydrophobic membrane ensures internal sterility prior to use. During priming, the FlowStop protection cap remains attached to the Male Luer-Lock adapter and allows air to pass but as soon as solution reaches the membrane, the flow stops automatically.



Disclaimer:

The above is indicated only for information. This document is not a specification. The values given are subject to the usual fluctuations. Please follow the instructions for use by the medication manufacturer and the actual SPC (Summary of Product Characteristics). Common hygiene measures and working methods are to be followed at all times.

Subject to alterations without prior notice!

REF	Product description	Packaging	GTIN	Filling/Storage volume
43.3030	L86 closed cap, 175 cm PVC (TOTM) tubing 3,0/4,1 mm (ID/OD), MLL adapter, FlowStop	1	04055698001293	7 ml/1 m (40°C)
		100	04055698001309	9 ml/1 m (40°C/2 bar)
		3600	-	
43.4975	L86 closed cap, 210 cm PVC (TOTM) tubing 3,0/4,1 mm (ID/OD), MLL adapter, FlowStop	1	04055698001736	7 ml/1 m (40°C)
		100	04055698001743	9 ml/1 m (40°C/2 bar)
		3600	-	

Product data sheet

UNSPSC: 42221609



Infusion - L86 - vented sets

L86 (closed cap, MLL)



REF	Product description	Packaging	GTIN	Filling/Storage volume
43.4349	L86 closed cap, 250 cm PVC (TOTM) tubing 3,0/4,1 mm (ID/OD), MLL adapter, FlowStop	1 100 3600	04055698001538 04055698001545 -	7 ml/1 m (40°C) 9 ml/1 m (40°C/2 bar)