

Light Inco Men

Absorption: LEVEL 2



Family:	Light Inco Men	Brand:	Egosan
Size:	One Size	Absorption Level:	Extra (4/10 drops)
Item Code:	XSU21B-14		
EAN-13 Code:	8011809015333	ITF-14 Code:	18011809015330
PRODUCT FEATURES			
Total weight (g/pcs, +/-5g):	22		
Overall dimensions finished product length x front width x back width (+/-20mm):	290x230x100		
Absorbent core dimensions (+/-20mm)			
Total length:	250		
Front side width:	180		
Back side width:	70		
Total absorption capacity ISO/DIS 11948-1 (g, +/-10%):	500		
PRODUCT SPECIFICATION			
Absorbent pad:	Pure ECF (Chlorine Free) cellulose and super absorbent polymer in sodium polyacrilate. ODOUR SECURE. Does not contain PHTHALATES		
Top sheet:	Soft hydrophilic non-woven, integrated with acquisition distribution layer		
Back sheet:	Breathable textile backsheet	(White)	
Individually wrapped:	Absent		
Glue:	Hypoallergenic, specific for medical devices		
Anti-leakage system	Soft elasticated leg cuffs		
Elastic systems	Latex free elastic threads (for leg cuffs)		
Adhesive Tapes:	Adhesive strip on the back of the product with relise paper		
Production batch:	Printed on the back		
LOGISTIC DATA			
Bag dimensions (mm - LxWxH):	230 x 90 x 180	N° Pieces per bag:	14
Bag gross weight (g):	320		
Box dimensions (mm - LxWxH):	580 x 240 x 365	N° bags per box:	12
Box gross weight (kg):	4,1		
N° of box per pallet:	21	Pallet gross weight (kg):	106,1
N° Layers per pallet :	3	Pallet protection:	wrapped with cellophane

Expire date: 5 years since date of production

Latex-free product

SANTEX S.p.A. - Via San Giovanni sul muro, 1 - 20121 Milano (MI) ITALY

Doc. n. 0 Rev n. 0 of 01/09/2023