



3-way stopcock- and Manifold extension lines

STERITEX® E25 (3W, SELF-LOCK)



REF	Product description	Packaging	GTIN	Filling/Storage volume
71.4507	30 cm PVC (TOTM) extension line 2,5/4,1 mm (ID/OD), 3W (ivory), SELF-LOCK, FlowStop	1	04055698034222	5 ml/1 m (40°C)
		80	04055698034239	5 ml/1 m (40°C/2 bar)
		5040	-	
71.4527	30 cm PVC (TOTM) extension line 2,5/4,1 mm (ID/OD), 3W (blue), SELF-LOCK, FlowStop	1	04055698034406	5 ml/1 m (40°C)
		80	04055698034413	5 ml/1 m (40°C/2 bar)
		5040	-	
71.4508	60 cm PVC (TOTM) extension line 2,5/4,1 mm (ID/OD), 3W (ivory), SELF-LOCK, FlowStop	1	04055698034246	5 ml/1 m (40°C)
		80	04055698034253	5 ml/1 m (40°C/2 bar)
		5040	-	
71.4528	60 cm PVC (TOTM) extension line 2,5/4,1 mm (ID/OD), 3W (blue), SELF-LOCK, FlowStop	1	04055698034420	5 ml/1 m (40°C)
		80	04055698034437	5 ml/1 m (40°C/2 bar)
		5040	-	
71.4509	90 cm PVC (TOTM) extension line 2,5/4,1 mm (ID/OD), 3W (ivory), SELF-LOCK, FlowStop	1	04055698034260	5 ml/1 m (40°C)
		80	04055698034277	5 ml/1 m (40°C/2 bar)
		5040	-	
71.4529	90 cm PVC (TOTM) extension line 2,5/4,1 mm (ID/OD), 3W (blue), SELF-LOCK, FlowStop	1	04055698034444	5 ml/1 m (40°C)
		80	04055698034451	5 ml/1 m (40°C/2 bar)
		5040	-	
71.4531	150 cm PVC (TOTM) extension line 2,5/4,1 mm (ID/OD), 3W (blue), SELF-LOCK, FlowStop	1	04055698034468	5 ml/1 m (40°C)
		80	04055698034482	5 ml/1 m (40°C/2 bar)
		5040	-	