


Description:	FORCEPS
Final dimensions:	13 cm
Composition:	Colourant, Reinforced Polypropylene.
Latex:	Free
PVC:	Free
Re use:	No
CMR:	Free
Endocrine Disruptors:	Free
Colour:	blue
Sterilization method:	Ethylene Oxide
Shelf life:	60 months
CFU/g	N.A.
Sterility Assur. Level_ SAL	≤10 ⁻⁶

STERILE EO






isento
Latex
free

isento
PVC
free



Nr and Revision Date:	31 -- 2021/08/24
CE Class:	Ila
CE mark:	CE1639
Compliance with:	MDD 93/42/EEC amended by 2007/47/EC
Applicable norms:	EN_ISO10993-1.09_AC.10 ; EN_ISO_10993-18.2009 ; EN_ISO_10993-18_2020 ; EN_ISO_10993-18_2020 ; EN_ISO_14971_2019 ; EN_ISO_15223-1.2021 ; EN_ISO_20417_2021 ; ISTA Standard 2A ; ICH Q1A (2003) ; ASTM D4332_01 ; ASTM D4169 ; EN_ISO_11607_1_2020 ; EN_ISO_11607-2.2020 ; EN_ISO_14971_2019_A11 ; ISO_2759.2014 ; ISO_3037_2013 ; EN_ISO_11135.2014_A1.2019 ; CEN_ISO_TS_16775_2021 ; EN_ISO_13485.2016 ; EN_ISO_13485.2016_AC18 ; EN_556-1.2001_AC.2006 ; ISO_11135.2014 ; EN_ISO_11607_1_2017 ; EN_ISO_11607-2.2017 ; ISO_TS_16775.2014 ; ASTM_F1980_16 ; ISO_14644-1.2015 ; ISO_14644-5_2004 ; ISO_14698_1_2003 ; EN_ISO_14971.2012 ; EN_ISO_15223-1.2016 ; EN_1041.2008_A1.2013 ;
GMDN:	62489 - Dressing/utility forceps, tweezers-like, single use

Level	Type	Dimension	Stamp	Labels
 1	1 Flexible peel-open pack	80x190mm		
 2	35 Folding box	78x193x360mm		labeled
 3	560 Carton box	786x386x321mm		labeled

Basic UDI	56056221206002JN
	EAN 1th level/nivel:05605622229539
	EAN 2th level/nivel:05605622429533
	EAN Carton Box/Caixa:05605622029535

Important Notes:	The device is intended to be used on a single patient during one single procedure to prevent risks from contamination, including possible infections, and/or functional damage.
-------------------------	---

Storage instructions	Store in a dry place at room temperature and protect from sunlight. The product has a good stability and therefore, normal conditions of medical devices storage shall apply.
-----------------------------	---

Waste instructions	No special requirements. The waste must be disposed according to local legal regulations for the disposable of medical waste.
---------------------------	---

Intended use	For grasping and to hold and stabilize tissues. For holding swabs for absorption during surgery or woundcare procedures or for skin disinfection procedures-
---------------------	--

Any modification of the product as presented or any use out of its intended use are not of the manufacturer responsibility.

ISSUED BY: Quality Management Department, 2022/09/15

CONFIDENTIAL: The information contained herein is confidential and belongs to Bastos Viegas, s.a. It is provided for specific purposes and in connection with established agreements. It is not permitted disclosure, reproduction, copying or distribution, in whole or in part, for other purposes or to any third parties.
