



Pants Classic

Size: X-Large Level: Extra



Family:	Pants Classic	Brand:	Egosan
Size:	X-Large	Absorption Level:	EXTRA (7/10 drops)
Item Code:	XPUX21B-14		
EAN-13 Code:	8011809016651	ITF-14 Code:	18011809016658
PRODUCT FEATURES			
Total weight (g/pcs, +/-5%):	95		
Overall dimensions finished product length x width x crotch (+/-20mm):	960x880x250		
Absorbent core dimensions (+/-20mm)			
Total length:	530		
Front side width:	190		
Back side width:	190		
Central crotch width:	105		
Waist size (cm):	110/165 (OVER 90kg)		
Total absorption capacity ISO/DIS 11948-1 (g, +/-10%):	1.650		
PRODUCT SPECIFICATION			
Absorbent pad:	Pure ECF (Chlorine Free) cellulose and super absorbent polymer in sodium polyacrilate. ODOUR SECURE. Does not contain PHTHALATE		
Top sheet:	Soft hydrophilic non-woven		
Back sheet:	Breathable textile backsheet cotton feel	White	
Elasticated pants chassis	Hypoallergenic elasticated waist band made in soft nonwoven, plus several elastic threads		
Anti-leakage system	Soft elasticated leg cuffs		
Acquisition distribution layer:	High Loft hydrophilic acquisition layer		
Glue:	Hypoallergenic, specific for medical devices		
Elastic systems	Latex free elastic threads (for crotch, leg cuffs and elasticated chassis; blue color in the rear belt)		
Side welding	Tear-away side seams for easier removal		
Disposal tape:	Z-Tape adhesive strip at the back for closing after use easy disposal		
Wetness indicator	Central fading away blu ink printing, with production batch indications		
LOGISTIC DATA			
Bag dimensions (mm - LxWxH):	260 x 160 x 260	N° Pieces per bag:	14
Gross weight bag (g):	1350		
Box dimensions (mm - LxWxH):	500 x 300 x 355	N° bags per box:	4
Gross weight Box (kg):	6,0	Box Features:	havana color, in-line printing
N° of box per pallet:	36	Gross weight Pallet (kg):	239,2
N° Layers per pallet :	6	Pallet protection:	Wrapped with cellophane

Expire date: 5 years since date of production

Latex-free product

SANTEX S.p.A. - Via San Giovanni sul muro, 1 - 20121 Milano (MI) ITALY

Doc. n. 0 Rev n. 1 of

07/10/2024