


# Oxygen Catheters

## Fixed Compress

SINGLE - USE



MEDICAL EQUIPMENT



### Area of application

Oxygen catheters are used for direct administration of oxygen via the nasopharyngeal route.

The products are for single use only.

Product range					
Connector	I.D./O.D. mm	Total Length mm	Tube Colour	CH	Reference
Oxygen Catheters with Fixed Compress (Shore A 76)					
Male	2.00/3.33	400	Blue	10	Q192010
Male	2.45/4.00	400	Blue	12	Q192012
Male	2.95/4.67	400	Blue	14	Q192014

Sterilization
Sterile, Sterilized by ethylene oxide, Shelf life: 5 years

Certifications	
CE-marking	CE 0123
Product classification	Class IIA (EU)
Quality system certified to	ISO 13485, Medical Device Single Audit Program (MDSAP)
European Medical Device Nomenclature (EMDN)	R03010203 OXYGEN THERAPY NASAL CANNULAS
Global Medical Device Nomenclature (GMDN)	12855

Instructions for use	
Included?	Yes

Material specifications	
Contains latex?	No
Contains phthalates?	No
Contains cotton?	No
Contains paper?	Packaging only
Tube and Connector	Polyvinyl chloride
Compress	Polyurethane foam

# Oxygen Catheters

## Fixed Compress

SINGLE - USE



MEDICAL EQUIPMENT

Packaging information									
REF	Primary pack./Single peelpack		Secondary pack./Shelfbox		Tertiary pack./Shipper			Gross weight kg	Pallet Qty. pcs
	Qty. pcs	Dimensions, mm (LxW)	Qty. pcs	Dimensions, mm (LxWxD)	Qty. pcs	Dimensions mm (LxWxD)	Volume m <sup>3</sup>		
Q192010	1	480x33.5	100	500x175x140	600	520x370x440	0.082	8.5	4,800
Q192012	1	480x33.5	100	500x175x140	600	520x370x440	0.082	9.5	4,800
Q192014	1	480x33.5	100	500x175x140	600	520x370x440	0.082	10.5	4,800

GTIN Codes			
REF	Primary package	Secondary package	Tertiary package
Q192010	6944932769878	16944932769875	26944932769872
Q192012	6936881154256	16936881154253	26936881154250
Q192014	6944932769878	16944932769875	26944932769872

Packaging material	
Primary package	Coated paper, Extruded polypropylene film
Secondary package	Corrugated fiberboard
Tertiary package	Double-wall corrugated fiberboard










# Oxygen Catheters

## Fixed Compress

SINGLE - USE



MEDICAL EQUIPMENT

Storage and handling	
	This way up
	Keep dry
	Fragile, handle with care
	Keep away from sunlight
	Do not resterilize
	Do not re-use
	Do not use if package is damaged and consult instructions for use
 -15°C 49°C	Temperature limit (-15° to 49°C)
	Consult instructions for use or consult electronic instructions for use

**SPONSORED IN AUSTRALIA BY**

Unoquip Pty Ltd  
PO Box 1117  
Windsor VIC 3181

1300 55 80 79  
[info@unoquip.au](mailto:info@unoquip.au)

**DISTRIBUTED BY**

**UNOQUIP GMBH**  
Schützengraben 20  
8200 Schaffhausen  
Switzerland

[info@unoquip.com](mailto:info@unoquip.com)  
[www.unoquip.com](http://www.unoquip.com)



Well Lead Medical Co., Ltd.  
C-4 Jinhui Industrial Estate  
Hualong, Panyu Guangzhou  
511434 P.R. China



Shanghai Intl. Holding  
Corp. GmbH  
Eiffestrasse 80  
20537 Hamburg  
Germany



QS Engineering AG  
Erlenstrasse 31  
CH - 4106 Therwil  
Switzerland

**UK  
Responsible  
Person  
(UKRP)**

SUNGO Certification  
Company Limited  
3rd Floor, 70 Gracechurch  
Street, London, EC3V 274HR  
England, United Kingdom